

Pro 2

The Most Advanced 3D Printer in Dentistry

Quick Start Guide

Rev1 - 110323

Congratulations!

You're about to unbox the most advanced dental 3D printer ever made. This guide contains everything you need to know to get started with Pro 2.

We're so excited to see how you use it to impact patient lives!

Training and Support

SprintRay University provides personalized virtual and in-person courses designed to maximize your printing experience. Discover more at sprintrayuniversity.com

If you have questions about your new products, please reach out to SprintRay:



Visit
Chat
Email

support.sprintray.com
sprintray.com/chat
support@sprintray.com



Sign Up account.sprintray.com

Sign Up with SprintRay

Your SprintRay account is used for everything from sending print jobs, ordering designs, and purchasing consumables.



You'll need an active account to send print files to your Pro 2.

Order Consumables

Visit the SprintRay store to order resin, tanks, and accessories for your Pro 2.

Visit store.sprintray.com

Order Appliance Designs

Submit scan files and patient information, receive 3D designs that are ready to print. Available treatments include occlusal guards, surgical guides, clear aligners, hybrid dentures, crowns, removable dentures, smile designs, veneers, bleaching trays, and more.

Unbox

Standard Equipment:

Pro 2 ships with a standard build platform and resin tank, optimized for maximum throughput. The build platform, resin tank, and Optical Panel are preinstalled.

Optional Equipment: (Sold Separately)

Optional accessories for the Pro 2, like the Arch Kit, can be fitted to help customize the printing experience for your office. Visit store.sprintray.com



Pro 2 3D Printer



Standard Build Platform
Preinstalled



Standard Resin Tank
and Lid
Preinstalled



Optical Panel
Preinstalled



Resin Wiper



Print Removal Tool



Support Snipper



Cleaning Cloth



Gloves



Documentation




Power Cable

Create a Workspace

The physical space your printer occupies plays a vital role in its success. Dust, heat, and direct sunlight can drastically impact your system's success or failure.

For best results, set Pro 2 up near your other SprintRay devices for optimal workflow efficiency. Turn the leveling feet on the bottom of the device to fine-tune the position of the printer and ensure that it sits level.

Recommended Operating Temperature



WARNING!

Do not place your printer near a mill or other equipment that produces dust. The inner components of the printer are highly sensitive, and excessive dust accumulation can impact print quality and accuracy.

Do


- Level surface
- Clean workspace
- Away from other equipment
- Away from windows and sunlight
- Consistent temperature 70-85 F
- Well ventilated area



Level Surface
Optimal Environment

Don't

- Uneven surface
- Dusty, dirty workspace
- Near other lab equipment
- Near windows or in sunlight
- Large temperature swings
- Poor ventilation




Uneven Surface
Warm Environment

Connecting Your SprintRay

Pro 2 is designed to be used with an active internet connection. Ordering designs, sending print jobs, and automatic handoff between devices require a SprintRay account with an active internet connection.


Advanced security protocols may prevent Pro 2 from connecting to the internet. Contact your network administrator for assistance.



Wireless Connection Great for most offices

A wireless connection is recommended for most offices that use normal networking and security protocols. Once you start up your Pro 2, you'll be prompted to connect to a wireless network by selecting your network name and entering the password.


i A wireless router is required to establish a wireless connection.



OR Wired Connection Great for offices with high security

A wired connection is recommended if your wireless network uses advanced security protocols or if your Wi-Fi signal is poor. Once you start up Pro 2, select the LAN option to connect to the internet.

i A network switch or router is required to establish a wired connection.



Network Security May Block Connection

Firewalls or antivirus software designed to secure your network may prevent your computer from finding the printer and connecting to SprintRay Cloud. If you need help getting connected, work with your network administrator and SprintRay support.

Hello, SprintRay


Pro 2 comes out of the box with the build arm lowered to the resin tank. The lid may also be locked to protect it in transit.

i When you turn Pro 2 on, the build platform will rise and the lid will unlock.

Plug In and Power Up

Use the provided power cable to plug your Pro 2 into an uninterrupted power supply or surge protector. To turn it on, press the power button.

When Pro 2 is on, press the power button to power off or restart. You can also press and hold the power button for 5 seconds to power it off.




Build Platform Lock

Pro 2 uses a magnetic lock to perfectly align and hold the build platform into place.



Unlocked
Remove to clean



Locked
Ready to print


Removal

To remove the build platform, turn the knob counterclockwise to unlock, then pull the platform toward you.

Installation


To install, slide the platform onto the arm, then turn clockwise to lock it in place.

Build Platform



Maintenance

Keep the platform clean and don't gouge it with sharp objects when removing prints. Do not drop.



Platform Heating

The build platform is equipped with a powerful heater. It may be warm to the touch.

Offline Printing


Use this USB-C port if you're having connectivity issues. Download print jobs onto a USB drive, then plug into Pro 2 and start printing.

Command Center

The ultrawide screen and powerful onboard computer allow you to create and manage print jobs, track the status of your consumables, and get step-by-step instructions for many workflows.

Sign In with SprintRay

To get the most out of your Pro 2, make sure to connect it to the internet and sign in with your SprintRay account.




Ergonomic Handle

Use the handle to open and close the lid to keep fingerprints and resin marks off the UV lid. Clean the handle with IPA 91% and a soft towel.


Resin Tank CONSUMABLE

Liquid resin is poured directly into this tank. Take note of the 'Max' and 'Min' lines - overfilling can cause damage to your printer. Resin tanks are consumable products. How long they last depends on the type of appliances printed, the resin used, and how they are maintained.



Tank Lid


This lid should be removed before each print job and reinstalled when you're finished printing for the day. This helps preserve the life of your resin.



Consumable Product

Your Pro 2 monitors the resin, usage, and health of your resin tank.

The onboard computer will notify you on the touchscreen when it's time to replace the tank.



Optical Panel REPLACEABLE

Located just beneath the resin tank, the Optical Panel is the heart of Pro 2.

Your printer monitors the usage and health of the Optical Panel.

i Learn more about when and how to replace the Optical Panel later in this guide.

Barcode Scanner

Pro 2 uses barcode scanning to track and log resin usage, expiration dates, and other important material information.

When you add or assign resin to a resin tank, the touchscreen will provide instructions for using the scanner.

Resin Tank CONSUMABLE

The resin tank holds liquid resin during the printing process. When you're not printing, cover the tank with the included lid.



Wear Gloves

Always wear gloves when handling uncured resin.



Do

- Use one tank for each resin type
- Keep away from light
- Hold by the handles
- Cover when not using
- Wear gloves while handling




Don't

- Change the resin type often
- Touch the film or underside
- Use sharp objects
- Puncture the film

Use the handles

The resin tank is designed with handles on either side. When you handle the resin tank, only hold it by the handles.




Pour Resin into the Tank

To activate a new resin tank, follow the instructions on the touchscreen to scan the bottle or select the resin type manually.

Pour resin directly from the bottle into the tank, filling until the resin is between the 'Min' and 'Max' lines.


i Do not overfill.



Check Tank Before Printing

Before each print job, use the included resin squeegee to stir the resin and check the bottom of the tank for cured resin or other debris. Mix until the resin is uniform in color.

i Mix until the resin is uniform in color.



Replace When Directed

The touchscreen on Pro 2 monitors resin usage and will notify you when it's time to replace a resin tank.

Failure to replace your resin tank in a timely manner may cause catastrophic damage to your printer that may not be covered under the warranty.

i For best results, keep an extra tank on hand and don't switch resin types.



Optical Panel REPLACEABLE

The Optical Panel is a long-term consumable item that is easy to replace in your office. You can track its health and usage via the printer touchscreen.

i Only replace the Optical Panel when directed by the printer touchscreen or customer support.

Keep it Clean

The Optical Panel should not require regular interaction, but it must remain clean at all times. Make sure there is no dust, resin, fingerprints, or other debris that might obstruct the path of light through the panel. Clean the top surface only - never touch the lens beneath.

Replaceable Product

The Optical Panel is a long-term consumable that must be replaced after it reaches a usage threshold that is monitored by the printer. To avoid potential downtime, keep an extra Optical Panel on hand. Purchase at store.sprinray.com

✓ Do



- Keep clean from dust and resin
- Follow instructions on the screen
- Wipe resin and IPA off immediately
- Hold by the front and sides

✗ Don't

- Touch the lens beneath
- Use abrasive cleaners
- Let IPA stand on the surface
- Use sharp objects on the panel

How to Clean

To clean, use IPA 91% and a soft cloth to wipe away any resin. Finish with glass cleaner to remove all streaks and dust.







Replacing Optical Panel

When the printer notifies you that it's time to change the Optical Panel, follow these instructions.




Wear Gloves

Always wear gloves when handling the Optical Panel.



The 3D Printing Workflow

This guide provides a general overview of the printing workflow, but the best way to learn is by doing. When you first boot up your printer, you'll be prompted to complete the first-time print, which will take you through each step in detail.



Add Resin and Stir

Install a resin tank that contains the material you'd like to print with. Make sure it is filled between the 'Min' and 'Max' lines in the tank.


Use the provided squeegee to mix the resin and make sure there is no hardened resin in the tank.

Check the Build Platform

Take off the build platform and run the scraper along the surface to remove any bits of uncured resin. Reinstall the platform and make sure it is locked into place.

Close the UV Lid

Keep the red hood closed when you're not directly accessing the build area.



Send Files to Printer


The best way to send print jobs to your Pro 2 is through SprintRay Cloud. You can build your own jobs with RayWare or directly print designs ordered through our Cloud Design service.

Visit account.sprinray.com to see your options. Once you've sent a job, check the printer touchscreen for further instructions.

Start the Print Job

Use the touchscreen to find or build the print job you'd like to run, then follow the onscreen instructions to complete the final checks.

Once you start the print job, the lid will lock and the platform will begin to heat the resin. You can follow the progress on the touchscreen.




Remove the Build Platform

When the print job is finished, the build platform will rise up to the resting position. Turn the knob counter-clockwise to unlock the build platform, then pull it toward you to remove.

Remove the Printed Part

If you're printing with supports, you should be able to remove the model from the build platform by hand. If you're printing something flat, like a dental model, use the provided scraper to remove the print.



Clean the Build Platform

Use the scraper to remove any hardened residue left on the build platform, then spray the surface with IPA and wipe it down with a soft cloth. Make sure there is no debris.

Reinstall the Build Platform

Once it's clean, slide the build platform over the build arm, then turn the knob clockwise to lock it into place.

Cover the Resin Tank

If you're finished printing for the day, secure the included resin tank cover over the tank to help preserve the life of the resin.

End-to-End Workflow

Pro 2 is designed to work best with SprintRay washing and curing devices for automatic handoff and processing profiles.

Visit sprinray.com to learn more about post processing devices.

